Table 1: Primary antibodies for immunohistochemistry.

Primary antibody	Dilution	Incubation time
Rabbit monoclonal anti-CD3, clone 2GV6,	Pre-diluted	16 min
Ventana Medical Systems, Tucson, AZ		
Rabbit monoclonal anti-CD4,clone SP35,	Pre-diluted	32 min
Ventana Medical Systems		
Rabbit monoclonal anti CD-8, clone SP57,	Pre-diluted	32min
Ventana Medical Systems		
Mouse monoclonal anti-CD20, clone L26,	Pre-diluted	16 min
Ventana Medical Systems		
Mouse monoclonal anti-CD68, clone KP-1,	Pre-diluted	12 min
Ventana,		
Mouse anti-DC-LAMP/CD208, clone 104G4,	1:50	32 min
#DDX0190, Dendritics, Dardilly, France		
Rabbit monoclonal anti-PD-L1, clone SP263,	Pre-diluted	16 min
#790-4905, Ventana Medical Systems		
Rat anti-PNAd, clone MECA-79, Biolegend, San	1:25	30 min
Diego, CA		
	Rabbit monoclonal anti-CD3, clone 2GV6, Ventana Medical Systems, Tucson, AZ Rabbit monoclonal anti-CD4,clone SP35, Ventana Medical Systems Rabbit monoclonal anti CD-8, clone SP57, Ventana Medical Systems Mouse monoclonal anti-CD20, clone L26, Ventana Medical Systems Mouse monoclonal anti-CD68, clone KP-1, Ventana, Mouse anti-DC-LAMP/CD208, clone 104G4, #DDX0190, Dendritics, Dardilly, France Rabbit monoclonal anti-PD-L1, clone SP263, #790-4905, Ventana Medical Systems Rat anti-PNAd, clone MECA-79, Biolegend, San	Rabbit monoclonal anti-CD3, clone 2GV6, Pre-diluted Ventana Medical Systems, Tucson, AZ Rabbit monoclonal anti-CD4, clone SP35, Pre-diluted Ventana Medical Systems Rabbit monoclonal anti CD-8, clone SP57, Pre-diluted Ventana Medical Systems Mouse monoclonal anti-CD20, clone L26, Pre-diluted Ventana Medical Systems Mouse monoclonal anti-CD68, clone KP-1, Pre-diluted Ventana, Mouse anti-DC-LAMP/CD208, clone 104G4, 1:50 #DDX0190, Dendritics, Dardilly, France Rabbit monoclonal anti-PD-L1, clone SP263, Pre-diluted #790-4905, Ventana Medical Systems Rat anti-PNAd, clone MECA-79, Biolegend, San 1:25